

# SAFETY WALKER

Model GB-510



# SAFETY WALKER Model GB-510

Helps seniors people to build self-reliance, and prevents them from being bedridden.

## Specifications

External dimensions	Approx. 630 (W) × 800 (L) × 890 to 1,210 (H) mm
Weight	Approx. 20 kg
Lifting method of seat	Gas spring
Lifting range of seat	Approx. 370 to 780 mm* stepless adjustment
Height adjustment of arm rest	Approx. 830 to 1,150 mm* stepless adjustment
Max. load weight	880 N (90 kgf) to the saddle 440 N (45 kgf) to each armrest (to vertically downward direction)

\*A little variety of height is expected because of the cushion.

Made in Japan



## Optional accessory Seat Cover Model GB510-001

### Specifications

External dimensions	Approx. 190 (W) × 350 (L) × 20 (H) mm
Weight	Approx. 200 g

Made in Japan



GB510-001

## Features and benefits

### ■ Allows patients to stand up and walk safely

With supports for the upper body and a seat for the lower body, Safety Walker allows user to move around safely without falling. The seat provides support if the user loses their balance or gets tired.

### ■ Helps patients get out of their wheelchairs, preventing them from being bedridden

Safety Walker can be used to increase the mobility of seniors who are often kept wheelchair-bound due to fall risk, as well as in the early-stage rehabilitation of hip fracture patients. If the user gets tired they can sit on the seat and rest, thus promoting independence by increasing patients' total exercise time and ambulation range.

### ■ Early-stage rehabilitation

Enables patients to begin rehabilitation right away after debilitating events, with upper- and lower-body supports to enable both standing practice and gait training.

### ■ Easy transfer

The seat height can be adjusted freely, in one motion. It can be lowered to the height of a bed or wheelchair, making it very easy to transfer patients to and from the Safety Walker. Also, because it's equipped with a gas spring mechanism to help patients stand up, even weaker patients can use it to stand up. The seat height can be adjusted using a quick and easy locking lever.

### ■ Removable seat

With a removable seat, Safety Walker is ideal for either standing practice or gait training, depending on the patient's functional ability.

### ■ Safe, wide casters

Safety Walker features wide casters to eliminate the risk of the casters getting stuck in gaps such as those in front of elevator doors, expanding the range in which users can move around safely.

- ① Upper body support with armrests and grips
- ② Rear fall guards
- ③ Stand-up support function with easy transfer from the bed or wheelchair
- ④ Arm rest lock levers
- ⑤ Wide casters with lock

- The colors in this brochure may not necessarily be the same as in the actual model due to printing conditions.
- Specifications and appearance are subject to change for improvement without notification.



OG Wellness Technologies Co., Ltd.

<http://www.og-wellness.com>

Head Office 1835-7 Miyoshi, Naka-ku, Okayama City, Okayama, 703-8261, JAPAN

TEL : +81-86-277-7181

FAX : +81-86-274-9072

Oku Factory 77 Mukoyama, Oku-cho, Setouchi City, Okayama, 701-4233, JAPAN